


# Connecting Roger to PC

1. If a camera is required, connect as normal and ensure it is selected under camera options, but no sound settings should be sent to it.
2. Locate the Roger Neckloop and Roger Multimedia hub.
3. Switch on your PC and the Roger devices.
4. Connect the Neckloop to the PC via the relevant USB to USB-c lead.
5. Connect the Multimedia hub via the headphone socket on your PC, using the 3.5mm jack to jack lead.



5. On your PC, **RIGHT** click on the speaker icon in the bottom right hand corner of your screen. Click 'Open Sound Settings'.
6. On the pop up window, ensure your PC's default option, '**Headphones**' is selected under Output. Ensure '**Microphone (Roger NeckLoop)**' is selected under Input.




# Connecting Roger to PC


The communication tool used in the following example is Microsoft Teams. Other apps may vary, but should be a similar process.

Although the sound settings on the PC are now configured, you may still need to configure the sound settings on Teams or your chosen app.

7. Click the 3 dots at the top. Select Settings. *On other apps this may be located elsewhere.*



8. Select Devices. Under Speaker, ensure that your PC's default option '**Headphones**' is selected. Under Microphone, ensure that '**Microphone (Roger NeckLoop)**' is selected.



9. The audio configuration should now be ready for use. On Teams, you can test this by doing a test call in the devices tab, under the Microphone options.